



Accesorios para Butirómetros

Gerber

Accesorios para Butirómetros



Code	Description	Packaging
LLG06202638	***no longer available*** ***replacemen	1 pz
LLG06224906	Rubber stopper double conical Ø 12/14x3	1 pz.
LLG06234037	Glass pin for milk butyrometer	1 pz.
LLG06262314	Glass beaker for butyrometers "butter"	1 pz.
LLG06264750	Weighing funnel, glass	1 pz.
LLG09112573	Rubber stopper GERBAL-M for volumetric m	1 pz.
LLG09112574	Key for GERBAL-M	1 pz.
LLG09112575	Rubber stoppers, large 30/17x22mm with h	1 pz.
LLG09112576	Rubber stoppers, small 15/10x13mm	1 pz.
LLG09112577	"Glass beaker for butyrometers" "cream" "	1 pz.
LLG09112578	"Glass beaker for butyrometers" "cheese" "	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

