



## Accesorios para garrafas serie 350

Accesorios para garrafas serie 350



| Code        | Description                               | Packaging |
|-------------|-------------------------------------------|-----------|
| LLG09011078 | Adapter, LDPE, blue with seal 6.303 821f  | 1 pz.     |
| LLG09011170 | Closure for stopcock opening with seal in | 1 pz.     |
| LLG09011181 | Stopcock, plastic, blue                   | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

