


**Tela metálica de acero inoxidable****RSG****Tela metálica de acero inoxidable**

Malla de acero inoxidable con bordes doblados. Luz de malla: 0,6 x 2,0 mm.



Code	Description	Packaging
LLG09150717	Wire ganges, stainless steel115x115mm,	1 pz.
LLG09150718	Wire ganges, stainless steel135x135mm,	1 pz.
LLG09150719	Wire ganges, stainless steel155x155mm,	1 pz.
LLG09150720	Wire ganges, stainless steel175x175mm,	1 pz.
LLG09150721	Wire ganges, stainless steel195x195mm,	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

