




**Tela metálica de acero inoxidable**

**RSG**

**Tela metálica de acero inoxidable**

Malla de acero inoxidable con bordes doblados. Luz de malla: 0,6 x 2,0 mm.



Code	Description	Packaging
LLG09150717	Wire ganges, stainless steel115x115mm,	1 pz.
LLG09150718	Wire ganges, stainless steel135x135mm,	1 pz.
LLG09150719	Wire ganges, stainless steel155x155mm,	1 pz.
LLG09150720	Wire ganges, stainless steel175x175mm,	1 pz.
LLG09150721	Wire ganges, stainless steel195x195mm,	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

