

**Probetas de vidrio, DURAN<sup>®</sup>****DURAN Group****Probetas de vidrio, DURAN<sup>®</sup>**

Con base. Reborde y superficie plana o reborde ligeramente esmerilado.

Esterilizables en autoclave.



Code	Description	Packaging
LLG09318021	Cylinders,DURAN,cap. 300 ml, 50x150 mm	1 pz
LLG09318034	Cylinders,DURAN,cap. 250 ml, 40x200 mm	1 pz
LLG09318036	Cylinders,DURAN,cap. 570 ml, 60x200 mm	1 pz
LLG09318046	Cylinders,DURAN,cap. 700 ml, 60x250 mm	1 pz
LLG09318052	Cylinders,DURAN,cap. 380 ml, 40x300 mm	1 pz
LLG09318053	Cylinders,DURAN,cap. 600 ml, 50x300 mm	1 pz
LLG09318068	Cylinders,DURAN,cap. 500 ml, 40x400 mm	1 pz
LLG09318074	Cylinders,DURAN,cap. 2000 ml, 80x400 mm	1 pz
LLG09318077	Cylinders,DURAN,cap. 1500 ml, 65x450 m	1 pz
LLG09318080	Cylinders,DURAN,cap. 1000 ml, 50x500 mm	1 pz
LLG09318107	Cylinders,DURAN,plane ground rim,cap.13	1 pz
LLG09318134	Cylinders,DURAN,plane ground rim,cap.25	1 pz
LLG09318136	Cylinders,DURAN,plane ground rim,cap.57	1 pz
LLG09318146	Cylinders,DURAN,plane ground rim,cap.70	1 pz
LLG09318168	Cylinders,DURAN,plane ground rim,cap.50	1 pz





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

