



Recipientes para muestras, PFA

Recipientes para muestras, PFA

Frascos de muestra de PFA con interior cónico y graduación elevada (subdivisiones de 5 ml).

Disponibles en dos tipos diferentes, dependiendo de la aplicación:

-Fondo exterior cónico

-Superficie plana en el fondo del frasco para mejorar la transferencia de calor (Recomendado para su uso con placas calientes)

Ambos tamaños de 50 ml caben en los estantes comunes de los automuestreadores.

El suministro se realiza sin tapón de rosca. Pedir el tapón de rosca (código 6.283 283) por separado.



Code	Description	Packaging
LLG06283278	Sample vessel 15ml, PFA plane, conica/1	1 pz.
LLG06283279	Sample vessel 25ml, PFA stand ring, c/1	1 pz.
LLG06283280	Sample vessel 25ml, PFA plane, conica/1	1 pz.
LLG06283281	Sample vessel 50ml, PFA stand ring, c/1	1 pz.
LLG06283282	Sample vessel 50ml, PFA plane, conica/1	1 pz.
LLG06283283	Screw cap, PFA for sample vessels 10//1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

