



Recipientes para muestras, PFA

Recipientes para muestras, PFA

Frascos de muestra de PFA con interior cónico y graduación elevada (subdivisiones de 5 ml). Disponibles en dos tipos diferentes, dependiendo de la aplicación:

- Fondo exterior cónico
- Superficie plana en el fondo del frasco para mejorar la transferencia de calor (Recomendado para su uso con placas calientes)

Ambos tamaños de 50 ml caben en los estantes comunes de los automuestreadores.

El suministro se realiza sin tapón de rosca. Pedir el tapón de rosca (código 6.283 283) por separado.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG06283278 | Sample vessel 15ml, PFA plane, conica/1 | 1 pz. |
| LLG06283279 | Sample vessel 25ml, PFA stand ring, c/1 | 1 pz. |
| LLG06283280 | Sample vessel 25ml, PFA plane, conica/1 | 1 pz. |
| LLG06283281 | Sample vessel 50ml, PFA stand ring, c/1 | 1 pz. |
| LLG06283282 | Sample vessel 50ml, PFA plane, conica/1 | 1 pz. |
| LLG06283283 | Screw cap, PFA for sample vessels 10//1 | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

