




Lupa de precisión de bolsillo easyPOCKET

Eschenbach

Lupa de precisión de bolsillo easyPOCKET
Color plata con iluminación: SMD-LED superplana,



Code	Description	Packaging
LLG09151147	Illuminated magnifying lens easyPOCKET3	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

