




Lupa de precisión de bolsillo easyPOCKET

Eschenbach

Lupa de precisión de bolsillo easyPOCKET

Color plata con iluminación: SMD-LED superplana,



Code	Description	Packaging
LLG09151147	Illuminated magnifying lens easyPOCKET3	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

