



Reactivos en solución para fotómetros de Lovibond<sup>®</sup>

Retsch

**Reactivos en solución para fotómetros de Lovibond<sup>®</sup>**

Solución de 15ml o botella de 100ml para rellenar el frasco de 15 ml.



Code	Description	Packaging
LLG06057943	Test solution DPD 1 15 ml flask UN 32/1	1 pz.
LLG06057944	Test solution DPD 3 flask 15 ml, red/1	1 pz.
LLG06058068	Buffer solution DPD 1 blue, flask 15m/1	1 pz.
LLG06203285	Test solution DPD 3 100ml flask/1	1 pz.
LLG07510822	Buffer solution DPD 1 blue, 100 ml fl/1	1 pz.
LLG07510823	Test solution DPD 1 green, 100 ml fla/1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

