



Ejes de agitación con roscas M6, M8 y M10 para cabezales de agitación

**Ejes de agitación con roscas M6, M8 y M10 para cabezales de agitación**  
Acero fino 1.4301. Para roscas M6, M8 y M10.



Code	Description	Packaging
LLG09156100	stirring rod for shaft with M6 thread15	1 pz
LLG09156101	stirring rod for shaft with M6 thread30	1 pz
LLG09156102	stirring rod for shaft with M6 thread35	1 pz
LLG09156103	stirring rod for shaft with M6 thread40	1 pz
LLG09156104	stirring rod for shaft with M6 thread45	1 pz
LLG09156105	stirring rod for shaft with M6 thread50	1 pz
LLG09156106	stirring rod for shaft with M6 thread60	1 pz
LLG09156107	stirring rod for shaft with M6 thread70	1 pz
LLG09156108	stirring rod for shaft with M6 thread80	1 pz
LLG09156129	stirring rod for shaft with M10 thread5	1 pz
LLG09156130	stirring rod for shaft with M10 thread7	1 pz
LLG09156131	stirring rod for shaft with M10 thread1	1 pz
LLG09156132	stirring rod for shaft with M10 thread1	1 pz
LLG09156133	Rod M 10 for stirrer heads550 x 12 mm,	1 pz
LLG09156134	Rods M 10 for stirrer heads750 x 12 mm,	1 pz
LLG09156135	Rods M 10 for stirrer heads800 x 12 mm,	1 pz
LLG09156136	Rods M 10 for stirrer heads1000 x 12 mm	1 pz
LLG09156137	Rods M 10 for stirrer heads1300 x 12 mm	1 pz
LLG09156138	Rods M 10 for stirrer heads1500 x 12 mm	1 pz





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

