



Varilla medidora del nivel de leche

Gerber

Varilla medidora del nivel de leche

latón niquelado, con válvula, esterilizable



Code	Description	Packaging
LLG06266491	Milk sample slave 1 ml d d	1 pz.
LLG06266492	Milk sample slave 2 ml brass nickel-pla	1 pz.
LLG06266493	Milk sample slave 4 ml brass nickel-pla	1 pz.
LLG06266494	Milk sample slave 5 ml brass nickel-pla	1 pz.
LLG06266495	Milk sample slave 10 ml brass nickel-pl	1 pz.
LLG06266496	Milk sample slave 20 ml brass nickel-pl	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

