



## Fascos de Woulff, DURAN<sup>®</sup>

### Fascos de Woulff, DURAN<sup>®</sup>

- Paredes gruesas, para trabajos con vacío
- Con recubrimiento de plástico como protección contra astillas e implosión
- Adaptador de vidrio con olivas roscadas de plástico (PP) Ø 11 mm
- Separador con válvula de ventilación
- Manómetro con 2 escalas (1.000 a 0 mbar, 760 a 0 mm Hg)



Code	Description	Packaging
LLG09305340	Woulff bottle 500ml, Erlenmeyer pattern	1 pz.
LLG09305341	Woulff bottles 1000ml, Erlenmeyer patte	1 pz.
LLG09305342	Woulff bottle 2000ml, Erlenmeyer patter	1 pz.
LLG09305343	Woulff bottles with pressure gange5000m	1 pz.
LLG09305344	Woulff bottles with pressure gauge10000	1 pz.
LLG09305345	Woulff bottles with pressure gange15.00	1 pz.
LLG09305346	Woulff bottles with pressure gange20.00	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

