



## Fascos de Woulff, DURAN<sup>®</sup>

### Fascos de Woulff, DURAN<sup>®</sup>

- Paredes gruesas, para trabajos con vacío
- Con recubrimiento de plástico como protección contra astillas e implosión
- Adaptador de vidrio con olivas roscadas de plástico (PP) Ø 11 mm
- Separador con válvula de ventilación
- Manómetro con 2 escalas (1.000 a 0 mbar, 760 a 0 mm Hg)



| Code        | Description                             | Packaging |
|-------------|---|-----------|
| LLG09305340 | Woulff bottle 500ml, Erlenmeyer pattern | 1 pz.     |
| LLG09305341 | Woulff bottles 1000ml, Erlenmeyer patte | 1 pz.     |
| LLG09305342 | Woulff bottle 2000ml, Erlenmeyer patter | 1 pz.     |
| LLG09305343 | Woulff bottles with pressure gange5000m | 1 pz.     |
| LLG09305344 | Woulff bottles with pressure gauge10000 | 1 pz.     |
| LLG09305345 | Woulff bottles with pressure gange15.00 | 1 pz.     |
| LLG09305346 | Woulff bottles with pressure gange20.00 | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

