



Adaptador de cámara C-Mount para la serie BA, AE & SMZ-171

Motic

Adaptador de cámara C-Mount para la serie BA, AE & SMZ-171



| Code | Description | Packaging |
|-------------|------------------------------|-----------|
| LLG06243486 | C-Mount camera adapter 0,65X | 1 pz |
| LLG07626193 | Adapter C-Mount 0,5x | 1 pz |



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

