

**Placas para el cultivo celular, tratadas, PS, estériles****Jet Bio-Filtration****Placas para el cultivo celular, tratadas, PS, estériles**

Estas placas de cultivo celular tratadas superficialmente de forma estándar presentan una evaporación reducida debido a unos rebordes de pocillo elevados y a unos aros uniformes en las tapas.

- Disponibles con 4-, 6-, 12-, 24-, 48-, 96- o 384-pocillos (tapa incluida)
- Pocillos con base plana
- Placa de 96 pocillos también disponible como placa de base redonda (6.267 171)
- Tapa de una sola posición que reduce los riesgos de contaminaciones cruzadas
- Los pocillos están etiquetados con un código alfanumérico para una identificación fácil
- Esterilizadas mediante radiación gamma
- No pirogénicas
- Material: Poliestireno



Code	Description	Packaging
LLG06267166	Tissue Culture plates 6-well, PS 10/100	100 uds
LLG06267167	Tissue Culture plates 12-well, PS 4/100	100 pz
LLG06267168	Tissue Culture plates 24-well, PS 2/100	100 uds
LLG06267169	Tissue Culture plates 48-well, PS 0/100	100 PZ
LLG06267170	Tissue Culture plates 96-well, PS 0/100	100 pz
LLG06267171	Tissue Culture plates 96-well, PS 0/100	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

