



Sonda para controladores de temperatura, TEMPAT<sup>®</sup>

**Sonda para controladores de temperatura, TEMPAT<sup>®</sup>**  
Todas las sondas con 2m de cable de conexión y enchufe.



Code	Description	Packaging
LLG07601585	Probe Pt 100, 3 mm Ø, 250 mm	1 pz
LLG09725320	Probe Pt 100, Ø 1,6 mm, 250 mmwith 2 m	1 pz
LLG09725321	Probe Pt 100, Ø 2,0 mm, 250 mmwith 2 m	1 pz
LLG09725322	Glass-probe Pt 100, Ø 2,5 mm250 mm	1 pz
LLG09725323	Glass-probe Pt 100, Ø 3,0 mm250 mm	1 pz
LLG09725394	Sensori per TEMPAT,sonda Fe-CuNi	1 pz
LLG09725395	Sensori per TEMPAT,sonda NiCr-Ni	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

