

**Tiras reactivas QUANTOFIX®****MACHERY NAGEL**

Las tiras reactivas QUANTOFIX® satisfacen todas las exigencias de los tests rápidos modernos. Las tiras reactivas son muy apreciadas por los usuarios profesionales para realizar análisis rápidos directamente en el lugar de toma de la muestra. Todos los tests QUANTOFIX se suministran listos para el uso. Las tiras reactivas han sido diseñadas para un solo uso, por lo que no necesitan mantenimiento ni calibración.



| Code     | Description                              | Packaging  |
|----------|--|------------|
| MNL91301 | QUANTOFIX Chromate box of 100 test stick | 100 strips |
| MNL91303 | QUANTOFIX Cobalt box of 100 test sticks  | 100 strips |
| MNL91304 | QUANTOFIX Copper box of 100 test sticks  | 100 strips |
| MNL91305 | QUANTOFIX Nickel box of 100 test sticks  | 100 strips |
| MNL91306 | QUANTOFIX Sulphite box of 100 test stick | 100 strips |
| MNL91307 | QUANTOFIX Aluminium box of 100 test stic | 100 strips |
| MNL91309 | QUANTOFIX Tin box of 100 test sticks 6x9 | 100 strips |
| MNL91310 | QUANTOFIX Zinc box of 100 test sticks 6  | 100 strips |
| MNL91311 | QUANTOFIX Nitrite box of 100 test sticks | 100 strips |
| MNL91312 | QUANTOFIX Peroxide 100 box of 100 test s | 100 strips |
| MNL91313 | QUANTOFIX Nitrate/Nitrite box of 100 tes | 100 strips |
| MNL91314 | QUANTOFIX Ascorbic acid box of 100 test  | 100 strips |
| MNL91315 | QUANTOFIX Ammonium box of 100 test stick | 100 strips |
| MNL91316 | QUANTOFIX Potassium box of 100 test stic | 100 strips |
| MNL91317 | QUANTOFIX Chlorine box of 100 test stick | 100 strips |
| MNL91318 | QUANTOFIX Cyanide box of 100 test sticks | 100 strips |
| MNL91319 | QUANTOFIX Peroxide 25 box of 100 test st | 100 strips |
| MNL91320 | QUANTOFIX Phosphate box of 100 test stic | 100 strips |
| MNL91321 | QUANTOFIX Chloride box of 100 test stick | 100 strips |
| MNL91322 | QUANTOFIX Nitrite 3000 box of 100 test s | 100 strips |



| Code     | Description                              | Packaging  |
|----------|--|------------|
| MNL91323 | QUANTOFIX Carbonate alkalinity box of 10 | 100 strips |
| MNL91324 | QUANTOFIX Calcium box of 60 test sticks  | 60 strips  |
| MNL91325 | QUANTOFIX Molybdate box of 100 test stic | 100 strips |
| MNL91326 | QUANTOFIX Multistick for Aquarium for de | 100 strips |
| MNL91327 | QUANTOFIX Multi-stick for aquarium owner | 25 strips  |
| MNL91328 | QUANTOFIX Formaldehyde pack of 100 test  | 100 strips |
| MNL91329 | QUANTOFIX Sulphate pack of 100 test stic | 100 strips |
| MNL91330 | QUANTOFIX Total iron 1000 Dip <(>&<)> Re | 100 strips |
| MNL91332 | QUANTOFIX Arsenic 50 pack of 100 test st | 100 strips |
| MNL91333 | QUANTOFIX Peroxide 1000 box of 100 test  | 100 strips |
| MNL91334 | QUANTOFIX Arsenic 10 for the determinati | 100 strips |
| MNL91335 | QUANTOFIX EDTA box of 100 test sticks 6x | 100 strips |
| MNL91336 | QUANTOFIX LubriCheck box of 100 test sti | 100 strips |
| MNL91337 | QUANTOFIX Quat box of 100 test sticks 6x | 100 strips |
| MNL91338 | QUANTOFIX Nitrit/pH for cooling lubrican | 100 strips |
| MNL91339 | QUANTOFIX Chlorine sensitive box of 100  | 100 strips |
| MNL91340 | QUANTOFIX Peracetic acid 50 box of 100 t | 100 strips |
| MNL91341 | QUANTOFIX Peracetic acid 500 box of 100  | 100 strips |
| MNL91342 | QUANTOFIX Peracetic acid 2000 box of 100 | 100 strips |
| MNL91343 | QUANTOFIX Glutaraldehyde box of 100 test | 100 strips |



| Code      | Description                              | Packaging   |
|-----------|--|-------------|
| MNL91344  | QUANTOFIX Total iron 100 box of 100 test | 100 strips  |
| MNL91345  | QUANTOFIX Arsenic Sensitive for the dete | 100 strips  |
| MNL91348  | QUANTOFIX Glucose box of 100 test sticks | 100 strips  |
| MNL91349  | QUANTOFIX active oxygen (potassium monop | 100 strips  |
| MNL91350  | QUANTOFIX Silver box of 100 test sticks  | 100 strips  |
| MNL91351  | QUANTOFIX Nitrate 100 apck of 100test    | 100 strips  |
| MNL91354  | Semi-quantitative test strips QUANTOFIX  | 100 test    |
| MNL91355  | QUANTOFIX Chlorine dioxide 15, 50 test s | 50 strips   |
| MNL91356  | Strisce QUANTOFIX Phosphate 10,strip 50  | 50 strips   |
| MNL91359  | Semi-quantitative test strips QUANTOFIX  | 25 pz       |
| MNL913918 | QUANTOFIX Nitrate set of 3 individually  | 3x50 strips |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

