

**Tiras reactivas QUANTOFIX®****MACHERY NAGEL**

Las tiras reactivas QUANTOFIX® satisfacen todas las exigencias de los tests rápidos modernos. Las tiras reactivas son muy apreciadas por los usuarios profesionales para realizar análisis rápidos directamente en el lugar de toma de la muestra. Todos los tests QUANTOFIX se suministran listos para el uso. Las tiras reactivas han sido diseñadas para un solo uso, por lo que no necesitan mantenimiento ni calibración.



Code	Description	Packaging
MNL91301	QUANTOFIX Chromate box of 100 test stick	100 strips
MNL91303	QUANTOFIX Cobalt box of 100 test sticks	100 strips
MNL91304	QUANTOFIX Copper box of 100 test sticks	100 strips
MNL91305	QUANTOFIX Nickel box of 100 test sticks	100 strips
MNL91306	QUANTOFIX Sulphite box of 100 test stick	100 strips
MNL91307	QUANTOFIX Aluminium box of 100 test stic	100 strips
MNL91309	QUANTOFIX Tin box of 100 test sticks 6x9	100 strips
MNL91310	QUANTOFIX Zinc box of 100 test sticks 6	100 strips
MNL91311	QUANTOFIX Nitrite box of 100 test sticks	100 strips
MNL91312	QUANTOFIX Peroxide 100 box of 100 test s	100 strips
MNL91313	QUANTOFIX Nitrate/Nitrite box of 100 tes	100 strips
MNL91314	QUANTOFIX Ascorbic acid box of 100 test	100 strips
MNL91315	QUANTOFIX Ammonium box of 100 test stick	100 strips
MNL91316	QUANTOFIX Potassium box of 100 test stic	100 strips
MNL91317	QUANTOFIX Chlorine box of 100 test stick	100 strips
MNL91318	QUANTOFIX Cyanide box of 100 test sticks	100 strips
MNL91319	QUANTOFIX Peroxide 25 box of 100 test st	100 strips
MNL91320	QUANTOFIX Phosphate box of 100 test stic	100 strips
MNL91321	QUANTOFIX Chloride box of 100 test stick	100 strips
MNL91322	QUANTOFIX Nitrite 3000 box of 100 test s	100 strips



Code	Description	Packaging
MNL91323	QUANTOFIX Carbonate alkalinity box of 10	100 strips
MNL91324	QUANTOFIX Calcium box of 60 test sticks	60 strips
MNL91325	QUANTOFIX Molybdate box of 100 test stic	100 strips
MNL91326	QUANTOFIX Multistick for Aquarium for de	100 strips
MNL91327	QUANTOFIX Multi-stick for aquarium owner	25 strips
MNL91328	QUANTOFIX Formaldehyde pack of 100 test	100 strips
MNL91329	QUANTOFIX Sulphate pack of 100 test stic	100 strips
MNL91330	QUANTOFIX Total iron 1000 Dip <(>&<)> Re	100 strips
MNL91332	QUANTOFIX Arsenic 50 pack of 100 test st	100 strips
MNL91333	QUANTOFIX Peroxide 1000 box of 100 test	100 strips
MNL91334	QUANTOFIX Arsenic 10 for the determinati	100 strips
MNL91335	QUANTOFIX EDTA box of 100 test sticks 6x	100 strips
MNL91336	QUANTOFIX LubriCheck box of 100 test sti	100 strips
MNL91337	QUANTOFIX Quat box of 100 test sticks 6x	100 strips
MNL91338	QUANTOFIX Nitrit/pH for cooling lubrican	100 strips
MNL91339	QUANTOFIX Chlorine sensitive box of 100	100 strips
MNL91340	QUANTOFIX Peracetic acid 50 box of 100 t	100 strips
MNL91341	QUANTOFIX Peracetic acid 500 box of 100	100 strips
MNL91342	QUANTOFIX Peracetic acid 2000 box of 100	100 strips
MNL91343	QUANTOFIX Glutaraldehyde box of 100 test	100 strips



Code	Description	Packaging
MNL91344	QUANTOFIX Total iron 100 box of 100 test	100 strips
MNL91345	QUANTOFIX Arsenic Sensitive for the dete	100 strips
MNL91348	QUANTOFIX Glucose box of 100 test sticks	100 strips
MNL91349	QUANTOFIX active oxygen (potassium monop	100 strips
MNL91350	QUANTOFIX Silver box of 100 test sticks	100 strips
MNL91351	QUANTOFIX Nitrate 100 apck of 100test	100 strips
MNL91354	Semi-quantitative test strips QUANTOFIX	100 test
MNL91355	QUANTOFIX Chlorine dioxide 15, 50 test s	50 strips
MNL91356	Strisce QUANTOFIX Phosphate 10,strip 50	50 strips
MNL91359	Semi-quantitative test strips QUANTOFIX	25 pz
MNL913918	QUANTOFIX Nitrate set of 3 individually	3x50 strips



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

