




Soporte para conos de sedimentación, PVC

Behr

Soporte para conos de sedimentación, PVC
Patrón de Imhoff, PVC gris sólido. Diseño estable.



Code	Description	Packaging
LLG09882101	Rack (holder) GSO 2for 2 Imhoff funnels	1 pz.
LLG09882102	Rack (holder) GSO 4for 4 Imhoff funnels	1 pz.
LLG09882103	Rack (holder) GSM 2for 2 Imhoff funnels	1 pz.
LLG09882104	Rack (holder) GSM 4for 4 Imhoff funnels	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

