



Vasos de molienda para molinos planetarios de bolas PULVERISETTE 5 / 6 classic line

Fritsch

**Vasos de molienda para molinos planetarios de bolas PULVERISETTE 5 / 6
*classic line***



Code	Description	Packaging
LLG09738315	Reducing adapter for beaker 80 ml	1 pz.
LLG09738512	Grinding jar,sintered carborundum 1,cap	1 pz.
LLG09738514	Grinding jar,Tungsten carbide,cap.80 ml	1 pz.
LLG09738515	Grinding beaker 80 ml, chrome steelwith	1 pz.
LLG09738518	Grinding jar,Zirconium oxide,cap.80 ml	1 pz.
LLG09738522	Grinding jar,sintered carborundum 1,cap	1 pz.
LLG09738524	Grinding jar,tungstene carbide,cap.250	1 pz.
LLG09738525	Grinding beaker, 250 ml, chrome steel h	1 pz.
LLG09738528	Grinding jar,Zirconium oxide,cap.250 ml	1 pz.
LLG09738532	Grinding jar,sintered carborundum 1,cap	1 pz.
LLG09738535	Grinding beaker 500 ml, chrome steel GE	1 pz.
LLG09738538	Grinding jar,Zirconium oxide,cap.500 ml	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

