



Tapas de porcelana para crisoles LLG, DIN 12904

LLG Labware

Tapas de porcelana para crisoles LLG, DIN 12904  
Vitrificadas.



Code	Description	Packaging
LLG09250930	LLG-Porcellaine lid D 30 DIN34 mm dia.	1 pz
LLG09250931	LLG-Porcelain lid D 35 DIN39 mm dia. fo	1 pz
LLG09250932	LLG-Porcelain lid D 40 DIN44 mm dia. fo	1 pz
LLG09250933	LLG-Porcelain lid D 45 DIN49 mm dia. fo	1 pz
LLG09250934	LLG-Porcellaine lid D 50 DIN54 mm dia.	1 pz
LLG09250935	LLG-Porcellaine lid D 60 DIN64 mm dia.	1 pz
LLG09250936	LLG-Porcelain lid D 70 DIN74 mm dia. fo	1 pz
LLG09250937	LLG-Porcellaine lid D 80 DIN84 mm dia.	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

