



Tapas de porcelana para crisoles LLG, DIN 12904

LLG Labware

Tapas de porcelana para crisoles LLG, DIN 12904
Vitrificadas.



Code	Description	Packaging
LLG09250930	LLG-Porcellaine lid D 30 DIN34 mm dia.	1 pz
LLG09250931	LLG-Porcelain lid D 35 DIN39 mm dia. fo	1 pz
LLG09250932	LLG-Porcelain lid D 40 DIN44 mm dia. fo	1 pz
LLG09250933	LLG-Porcelain lid D 45 DIN49 mm dia. fo	1 pz
LLG09250934	LLG-Porcellaine lid D 50 DIN54 mm dia.	1 pz
LLG09250935	LLG-Porcellaine lid D 60 DIN64 mm dia.	1 pz
LLG09250936	LLG-Porcelain lid D 70 DIN74 mm dia. fo	1 pz
LLG09250937	LLG-Porcellaine lid D 80 DIN84 mm dia.	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

