



Accesorios para molino discontinuo Tube Mill control / Tube Mill 100 control

Accesorios para molino discontinuo Tube Mill control / Tube Mill 100 control



Code	Description	Packaging
LLG04652716	Disposable grinding chamber MT 40.10/10	10 pz.
LLG04653944	Abrasion set A-MMT 40.100 for MMT 40.1	1 pz.
LLG06260495	Grinding beaker MT 40 pack of 100/100	100 pz.
LLG06285802	Disposable grinding jar 100ml MT 100.10,	10 pz.
LLG06310247	Disposable grinding jar 100ml MT 100.10,	10 pz.
LLG07639994	Multiple grinding chamber MMT 40.1	1 pz.
LLG07644728	Disposable grinding beaker MT 100.50 100	50 pz.
LLG09737791	Disposable grinding chamber MT 40pac/10	10 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

