



Accesorios para agitadores suspendidos Hei-TORQUE

Heidolph

Accesorios para agitadores suspendidos Hei-TORQUE



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG06228500 | Stand S2 XXL | 1 pz. |
| LLG07047043 | Flexible shaft including chuck, 8mm | 1 pz. |
| LLG07621811 | Telescopic stand | 1 pz. |
| LLG09816508 | Chuck 8,5mm span | 1 pz. |
| LLG09816562 | Support clamp for holder rod 13 mm | 1 pz. |
| LLG09816564 | Stand S2 | 1 pz. |
| LLG09816566 | Stirrer shaft guard, plexiglass, height | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

