



## Sistema de filtros compatible para máscaras completas y semimáscaras 3M

### Sistema de filtros compatible para máscaras completas y semimáscaras 3M™

Todos los filtros se adaptan a la conexión mediante clic a la serie de máscaras 6000 , 6500 , y 7500 y 7000. El filtro adecuado para cada aplicación. No importa cuál sea el área de trabajo, utilizando mascarillas de respiración 3M™ , el amplio programa de filtros 3M™ ofrece una seguridad óptima y proporciona posibilidades de aplicación versátiles. La gama de filtros compatibles 3M™ incluye filtros combinables e intercambiables contra gases, vapores y/o partículas y también filtros combinados.



Code	Description	Packaging
LLG06267784	Particle insertion filter P1R against/2	2 pz.
LLG06267785	Particle insertion filter P2R against/2	2 pz.
LLG06267786	filter P3R PACK 2	2 pz.
LLG06269584	Cartridge against Ammonia K1 pack of /2	2 pz.
LLG06269585	Gas and combined filters A1HgP3R for /2	2 pz.
LLG06269586	Filter for masks 6000/7000/7500 A1+fo/2	2 pz.
LLG06269587	Filter Lid for combination of gas filter	2 pz.
LLG06269588	Cartridge A2 against Organic Vapors	2 pz.
LLG06269589	Cartridge ABE1 against Organic Inorganic	2 pz.
LLG06269590	Cartridge A1 against Organic Vapors	2 pz.
LLG06270183	Gas- and combination filter AXP3 for /2	2 pz
LLG07072310	FILTRI PERMASCHERE SEMI E PIENO FACCIALI	2 pz
LLG07620305	Cartridge ABEK1 against Organic Inorgan	8 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

