



Mechero Bunsen Labogaz[®] 470

Mechero Bunsen Labogaz[®] 470

Para usarse con cartuchos de gas CV 300 Plus (código 6.243 905) o CV 470 Plus (código 9.018 512).

Peso: 380 g sin cartucho

Altura: 190 mm con cartucho CV 300 Plus

240 mm con cartucho CV 470 Plus

Consumo: 55 g / h



Code	Description	Packaging
LLG09018511	Bunsen burner,for use with gas cartridg	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

