



Accesorios para espectrofotómetro LLG-uniSPEC 2 / LLG-uniSPEC 4

Accesorios para espectrofotómetro LLG-uniSPEC 2 / LLG-uniSPEC 4



Code	Description	Packaging
LLG06263615	Tungsten lamp for uniSPEC 2 and uniSP/1	1 pz.
LLG06263616	Deuterium lamp for uniSPEC 2 and uniS/1	1 pz.
LLG06263617	Adjustable cell holder for 1 micro-cu/1	1 pz.
LLG06266530	Cell holder for 4 cuvettes 50mm for u/1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

