

**Gas safety tubing, rubber**

Flexible connecting line of a rubber tubing with a contracted metal spiral and double-sided end sockets according to DVGW 5501(P). Kink resistant. For connecting between safety gas connection fitting according to DIN 3383-4 and gas burners for laboratories according to DIN 30665 section 1. For all gases according to DVGW worksheet G 260, at pressures up to 100 mbar and temperatures up to 70°C. Installations according to DVGW-TRGI 2008 and DVGW worksheet G 621. DVGW Reg.-No. DG-4607AU2366. Sold only in specified length



Code	Description	Packaging
LLG09018935	Gas safety tubing, rubber, Length 500 mm	1 pz
LLG09018936	Gas safety tubing, rubber, Length 600 mm	1 pz
LLG09018937	Gas safety tubing, rubber, Length 750 mm	1 pz
LLG09018940	Gas safety tubing, rubber, Length 1000 mm	1 pz
LLG09018942	Gas safety tubing, rubber, Length 1250 mm	1 pz
LLG09018945	Gas safety tubing, rubber, Length 1500 mm	1 pz
LLG09018950	Gas safety tubing, rubber, Length 2000 mm	1 pz
LLG09018960	Gas safety tubing, rubber, Length 3000 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

