



Rotores para Concentrator plus

Eppendorf

Rotores para Concentrator plus™



Code	Description	Packaging
LLG09776465	Fixed angle rotor F-45-48-11for 48 tube	1 pz.
LLG09776490	Rotor for concentrator 5301for 8 tubes	1 pz.
LLG09776496	Rotor A-2-VC, 2 x 96x0.2 ml	1 pz.
LLG09776624	Rotor F-35-6-30for 6x15ml Falcon tubes	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

