

**Anillo de vertido para tapón de rosca Premium, TpCh260****DURAN Group****Anillo de vertido para tapón de rosca Premium, TpCh260**

Anillo de vertido incoloro TpCh260, altura 4mm. Resistente a temperaturas de -196°C a +260°C.



Code	Description	Packaging
LLG06803677	Premium-pouring ring GL 45height 4 mm,	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

