



Tapones para frascos de cuello ancho, LDPE

Kautex Textron GmbH & Co.KG

Tapones para frascos de cuello ancho, LDPE



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06291538 | Screw cap 40mm LDPE /wide mouth bottle | 1 pezzo |
| LLG06291539 | Screw cap 50 mm, LDPE for wide mouth bot | |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

