



Cestos para autoclaves verticales digitales, serie AES

Cestos para autoclaves verticales digitales, serie AES



Code	Description	Packaging
LLG06261937	Wire basket CV-75, stainless steel Ø /1	1 pz.
LLG06268955	Wire basket CV-12, stainless steel Ø /1	1 pz
LLG06270951	Wire basket for AES-28stainless steel,	1 pz.
LLG09851876	Solit basket for AE-AES-28stainless ste	1 pz
LLG09851877	Solid basket for AE-AES-75/110stainless	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

