



Accesorios para registrador de datos PDF EBI 300 / EBI 310

Accesorios para registrador de datos PDF EBI 300 / EBI 310



Code	Description	Packaging
LLG06263842	Humidistat TPH 500humidity: 0...100%0..	1 pz.
LLG06281569	Stainless Steel-Sinter-Filter AH 300/1	1 pz.
LLG09908017	Wall fastener EBI 300 WM2for EBI300 and	1 pz.
LLG09908027	External NTC probe TPC 300for EBI 300	1 pz.
LLG09908099	Temperature probe TPX 220Pt1000, -200..	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

