



Escurreidores, alambre recubierto de PE

Escurreidores, alambre recubierto de PE

Montaje en pared o sobre mesa. Con canal de goteo y patas con topes de plástico.



Code	Description	Packaging
LLG09003305	Draining racks,PE-coated wirefor 60 tes	1 pz.
LLG09003320	Draining racks,PE-coated wirefor 24 tes	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

