




## Soporte para portaobjetos, PE

### Soporte para portaobjetos, PE

Para portaobjetos de 76 x 26 mm. Para almacenar, limpiar o secar los portaobjetos. Diseño robusto con superficie de goma antideslizante. Resistente a temperaturas de -50 a +80 °C.



Code	Description	Packaging
LLG06254573	Slide holder 153x102x21mm PE, for 15 /1	1 pz.
LLG06254574	Slide holder 298x105x21mm PE, for 40 /1	1 pz.
LLG06254575	Slide holder 299x178x21mm PE, for 78 /1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

