




## Soporte para portaobjetos, PE

### Soporte para portaobjetos, PE

Para portaobjetos de 76 x 26 mm. Para almacenar, limpiar o secar los portaobjetos. Diseño robusto con superficie de goma antideslizante. Resistente a temperaturas de -50 a +80 °C.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG06254573 | Slide holder 153x102x21mm PE, for 15 /1 | 1 pz.     |
| LLG06254574 | Slide holder 298x105x21mm PE, for 40 /1 | 1 pz.     |
| LLG06254575 | Slide holder 299x178x21mm PE, for 78 /1 | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

