




Soporte para portaobjetos, PE

Soporte para portaobjetos, PE

Para portaobjetos de 76 x 26 mm. Para almacenar, limpiar o secar los portaobjetos. Diseño robusto con superficie de goma antideslizante. Resistente a temperaturas de -50 a +80 °C.



Code	Description	Packaging
LLG06254573	Slide holder 153x102x21mm PE, for 15 /1	1 pz.
LLG06254574	Slide holder 298x105x21mm PE, for 40 /1	1 pz.
LLG06254575	Slide holder 299x178x21mm PE, for 78 /1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

