



Recipientes Dewar, forma cilíndrica, para CO<sub>2</sub> y LN<sub>2</sub>.

KGW

**Recipientes Dewar, forma cilíndrica, para CO<sub>2</sub> y LN<sub>2</sub>.**

DURAN®. DIN 12492. Forma cilíndrica. Vidrio de borosilicato 3.3. Con revestimiento protector azul metálico o con estructura de aluminio. Gran abertura. Con tapa aislada. Asas de transporte como se indica. Para enfriar, almacenar y transportar pequeñas muestras de material sensible a la temperatura, por ejemplo, con hielo seco (CO<sub>2</sub>). Al almacenar nitrógeno líquido (LN<sub>2</sub>) es preciso un agujero en el centro de la tapa (diámetro aprox. 2mm).



Code	Description	Packaging
LLG09032131	Large insulating vessels,with cover and	1 pz.
LLG09032132	Large insulating vessels,with cover and	1 pz.
LLG09032133	Large insulating vessels,with cover and	1 pz.
LLG09032134	Large insulating vessels,with cover and	1 pz.
LLG09032742	Dewar flasks 10I31 CAL, bowl form	1 pz.
LLG09032743	Dewar flasks 21I33CAL, bowl form	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

