



Recipientes Dewar, forma cilíndrica, para CO₂ y LN₂.

KGW

Recipientes Dewar, forma cilíndrica, para CO₂ y LN₂.

DURAN®. DIN 12492. Forma cilíndrica. Vidrio de borosilicato 3.3. Con revestimiento protector azul metálico o con estructura de aluminio. Gran abertura. Con tapa aislada. Asas de transporte como se indica. Para enfriar, almacenar y transportar pequeñas muestras de material sensible a la temperatura, por ejemplo, con hielo seco (CO₂). Al almacenar nitrógeno líquido (LN₂) es preciso un agujero en el centro de la tapa (diámetro aprox. 2mm).



Code	Description	Packaging
LLG09032131	Large insulating vessels,with cover and	1 pz.
LLG09032132	Large insulating vessels,with cover and	1 pz.
LLG09032133	Large insulating vessels,with cover and	1 pz.
LLG09032134	Large insulating vessels,with cover and	1 pz.
LLG09032742	Dewar flasks 10I31 CAL, bowl form	1 pz.
LLG09032743	Dewar flasks 21I33CAL, bowl form	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

