



Tapones, serie 303, HDPE

Kautex

Tapones, serie 303, HDPE
Marrones, tapones de rosca.



Code	Description	Packaging
LLG09073778	Screw cap LDPE, brown for 50 and 100ml w	1 pz.
LLG09073779	LDPE screw-type cap brown for WH bottle	1 pz.
LLG09073781	Bell cap PP black ring red	1 pz.
LLG09073782	Screw caps, brown, LDPE for 1000 ml wide	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

