

**Mono para sala limpia ASPURE, cremallera lateral, poliéster****Mono para sala limpia ASPURE, cremallera lateral, poliéster**

El mono está disponible en cuatro versiones de color diferentes, sin capucha ni máscara. La caída del cabello, el polvo y otras partículas se evita con mangas dobles. En la zona de las mangas, se incorporan las bandas de goma de la cintura y del dobladillo inferior.

- Cremallera lateral oblicua
- Protección antiestática
- Con hilo conductor
- Adecuado para la clase de sala limpia ISO 5
- Resistencia de la superficie: 106 - 108 @?

Para el uso de la sala limpia, esta prenda debe ser lavada de acuerdo a la clase de sala limpia requerida.



Code	Description	Packaging
LLG04659040	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659041	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659042	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659043	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659044	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659045	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659046	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659047	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659048	ASPURE Overall for Cleanroom, lateral zi	1 pz.
LLG04659049	ASPURE Overall for Cleanroom, lateral zi	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

