

**Cierres para la serie 310 PETG y PVC****Kautex****Cierres para la serie 310 PETG y PVC**

Tapas de rosca negras, adecuadas para contenedores de boca ancha de PVC (transparente o ámbar) y PETG.



Code	Description	Packaging
LLG09073070	Screw cap, 23 mm	1 pz.
LLG09073075	Screw cap, 69 mm	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

