



## Tarros con tapón de rosca Nalgene, PP

### Tarros con tapón de rosca Nalgene™, PP

Transparentes. Excelente resistencia química. Autoclavables - deshacer la tapa antes del autoclave.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09033230 | Jars,PP,with screw cap,cap. 15mlpack/12 | 12 pz.    |
| LLG09033232 | Jars,PP,with screw cap,cap. 60 mlpac/12 | 12 pz.    |
| LLG09033233 | Jars,PP,with screw cap,cap. 125 mlpa/12 | 12 pz.    |
| LLG09033234 | Jars,PP,with screw cap,cap. 250 mlpac/6 | 6 pz.     |
| LLG09033235 | Jars,PP,with screw cap,cap. 500 mlpac/6 | 6 pz.     |
| LLG09033236 | Jars,PP,with screw cap,cap. 950 mlpac/6 | 6 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

