



Tarros con tapón de rosca Nalgene, PP

Tarros con tapón de rosca Nalgene™, PP

Transparentes. Excelente resistencia química. Autoclavables - deshacer la tapa antes del autoclave.



Code	Description	Packaging
LLG09033230	Jars,PP,with screw cap,cap. 15mlpac/12	12 pz.
LLG09033232	Jars,PP,with screw cap,cap. 60 mlpac/12	12 pz.
LLG09033233	Jars,PP,with screw cap,cap. 125 mlpa/12	12 pz.
LLG09033234	Jars,PP,with screw cap,cap. 250 mlpac/6	6 pz.
LLG09033235	Jars,PP,with screw cap,cap. 500 mlpac/6	6 pz.
LLG09033236	Jars,PP,with screw cap,cap. 950 mlpac/6	6 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

