

**Jars with screw cap Nalgene, PP**

Thermo Scientific

Jars with screw cap Nalgeneä, PP

Transparent. Excellent chemical resistance. Autoclavable - undo cap before autoclaving.



Code	Description	Packaging
LLG09033150	Jars with screw cap Nalgene, PP, Nominal capacity 15 ml	
LLG09033151	Jars with screw cap Nalgene, PP, Nominal capacity 30 ml	
LLG09033152	Jars with screw cap Nalgene, PP, Nominal capacity 60 ml	
LLG09033153	Jars with screw cap Nalgene, PP, Nominal capacity 125 ml	
LLG09033154	Jars with screw cap Nalgene, PP, Nominal capacity 250 ml	
LLG09033155	Jars,PP,with screw cap,cap. 500 ml/1	1 pz
LLG09033156	Jars with screw cap Nalgene, PP, Nominal capacity 1000 ml	
LLG09033230	Jars with screw cap Nalgene, PP, Nominal capacity 15 ml	12 pz
LLG09033231	Jars with screw cap Nalgene, PP, Nominal capacity 30 ml	12 pz.
LLG09033232	Jars with screw cap Nalgene, PP, Nominal capacity 60 ml	12 pz
LLG09033233	Jars with screw cap Nalgene, PP, Nominal capacity 125 ml	12 pz
LLG09033234	Jars with screw cap Nalgene, PP, Nominal capacity 250 ml	6 pz
LLG09033235	Jars with screw cap Nalgene, PP, Nominal capacity 500 ml	6 pz
LLG09033236	Jars with screw cap Nalgene, PP, Nominal capacity 1000 ml	6 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

