



Fuente de luz fría, accesorios, SCHOTT

Schott Fiber Optics

Fuente de luz fría, accesorios, SCHOTT



Code	Description	Packaging
LLG06270979	Foot switsch for VisiLED-Controller M/1	1 pz.
LLG07655134	Mounting clamp for cast base 158.340 /1	1 pz.
LLG09705196	Halogen lamp type EFR 15V/150W, for K/1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

