



Sondas de temperatura, Pt100 para ama-digit ad 20 th

Sondas de temperatura, Pt100 para ama-digit ad 20 th

Clase A. Sonda revestida de acero inoxidable, aprox. 105 x 3mm L x diám.

Con mango y aprox. 1m de cable.



Code	Description	Packaging
LLG09243111	Immersion sensor Pt-100	1 pz.
LLG09243112	Penetration sensor Pt-100E910652	1 pz.
LLG09243113	Surface sensor Pt-100	1 pz.
LLG09243114	Air sensor Pt-100	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

