



Viscosímetros con boquilla fija

Viscosímetros con boquilla fija

Con boquilla fija a base de acero inoxidable. La carcasa es de aluminio anodizado de titanio.



Code	Description	Packaging
LLG06091529	Flow cups DE 10, aluminium 3 mm flow no	1 pz.
LLG06226824	Dipping flow cups DE 10, DIN 53211 with	1 pz.
LLG06301841	Dipping flow cups DE 10, aluminum with	1 pz.
LLG06302759	Viscosity cup,DIN 53211,Aluminium,with	1 pz.
LLG09149862	Viscosity cup,DIN 53211,Aluminium,with	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

