



**Cestos de alambre, alambre/nylon**

**Cestos de alambre, alambre/nylon**

Blancos. Alambre resistente revestido de nylon.



| Code        | Description                             | Packaging |
|-------------|---|-----------|
| LLG09033341 | Wire baskets,nylon-coated,100x100x100 m | 1 pz.     |
| LLG09033342 | Wire baskets,nylon-coated,120x120x120 m | 1 pz.     |
| LLG09033344 | Wire baskets,nylon-coated,140x140x140 m | 1 pz.     |
| LLG09033346 | Wire baskets,nylon-coated,160x160x160 m | 1 pz.     |
| LLG09033350 | Wire baskets,nylon-coated,200x200x200 m | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

